

HyPro medical



A hydrogen peroxide solution to clean and disinfect in places where hygiene is vital.

Analisis commercializes and develops Hymetec products and equipment which act effectively on any kind of infection with the aim to obtain the best possible hygiene for everyone.

HyPro medical, H₂O₂ 3% solution, to efficiently clean and disinfect surfaces of non-invasive medical equipment and devices.

Technical information

Use	Ready-to use solution
Active substance	Hydrogen peroxide 3%
Tensioactive agent	Anionic
Color/form	Colorless-limpid
Perfumes	Green apple
GMDN	45058
Shelf Life	18 months unopened 2 months after opening
Stability	Stable from 5°C à 30°C
Relative density	1
Storage conditions	In original packaging keep in ventilated place. Keep away from direct sunlight. Keep bottles closed.
Packaging	Spray 1 L (boites de 6) 1 L foam (boxes of 6) Can 5 L

Efficiency

Tests were carried out according to Standards EN 13727 + A2 (bactericidal activity: *S. aureus*, *P. aeruginosa*, *E. hirae*) and EN 14476 (virucidal activity: polio virus 1).
HyPro medical eliminates more than 99.9% of the total flora.

Applications

HyPro medical is applicable in various fields such as :

- Hospitals
- Medical centers
- Ambulances
- Laboratories

Precautions of use

- When handling, please wear suitable protective clothes (safety glasses, chemical gloves and safety mask).
- Check the expiration date of the product.
- Please consult the safety data sheet (available on our website: www.analis.be/hymetec).
- Handle the product away from any heat.
- Do not reuse the empty bottle

Benefits

- ✓ No residue.
- ✓ No corrosion.
- ✓ Efficient on bacteria and viruses.
- ✓ Cleaning and disinfection of surfaces of non-invasive medical equipment and devices.
- ✓ Ready-to-use solution.
- ✓ Eco-friendly product.
- ✓ Does not contain quaternary ammonium, phenol, chlorine, aldehyde, alcohol and VOC.